

## OVERVIEW

Waverly's new C100 minifuge is the most compact, quiet and elegantly designed mini microcentrifuge in the world! A must-have personal centrifuge for all laboratory workers to provide convenient quick spin downs right at the bench.

The Waverly C100 Mini Centrifuges come complete with 2 rotors: 6 x 1.5/2.0ml microcentrifuge tubes, and 2 x 0.2ml 8-strip PCR tubes. Adapters are also included for 0.2ml and 0.5ml tubes. C100 is very simple to use, just close the lid and the centrifuge will reach top speed within seconds. When finished, open the lid and the rotor will automatically brake to a halt. This design maximizes time efficiency for all those lab applications that just need quick spin down separation.



## PRODUCT FEATURES

- Ideal for microfiltration and quick spin downs of microtubes tubes and PCR tubes
- Includes 2 rotors for 6 x 1.5/2.0ml microcentrifuge tubes and 2 x 8-strip 0.2ml PCR tubes
- Very quiet operation with noise level only  $\leq 35$  dB
- 5,400 RPM / 2,000 x g
- Fast Acceleration and Braking within seconds of opening/closing lid
- Easily change rotors without tools
- Elegant, compact, ergonomic design; fits in the palm of your hand
- 2-Year Warranty

## TECHNICAL SPECIFICATIONS

Max Speed	5,400 RPM
Max RCF	2,000 x g
Rotor Capacity	6 x 0.2-2.0ml tubes or 2 x 8-strip PCR tubes
Run Time	Continuous operation
Drive Motor	DC Motor
Power	AC 110-240V / 50-60Hz
Noise Level	$\leq 35$ db
Dimensions	6.9L x 5W x 4H inches
Weight	2.5 lbs.



## ORDERING INFORMATION

**C100**

Waverly C100 Mini Centrifuge, Clear Lid, 120V

## WHAT'S INCLUDED

- Waverly C100 microcentrifuge, 120V
- Rotor for 6 x 1.5/2.0ml microcentrifuge tubes
- Rotor for 16 x 0.2ml tubes (or 2 x 8-strip PCR tubes)
- Adapters for 0.5ml tubes