



Z 100 A

CLINICAL LAB CENTRIFUGE FOR 5ML, 7ML, 10ML AND 15ML TUBES

Labnet
Labnet International, Inc.




- Ideal for doctor's offices and small clinics
- Useful for blood, urine, and cell separations
- Supplied with 6 x 15ml rotor
- Variable speed from 1,000-4,000 rpm (1,900 x g)
- Capable of centrifuging all common blood draw tubes
 - 15ml (16 x 124mm)
 - 10ml (16 x 100mm)
 - 7ml (13 x 100mm)
 - 5ml (13 x 75mm)

SPECIFICATIONS

Speed control	Variable
Max. speed	4,000 rpm
Max. g - force	1,900 x g
Timer	1-30 min or continuous
Rotor angle	45°
Dimensions (WxDxH)	10.5 x 13.4 x 9.75 in 25 x 32 x 23 cm
Weight	12 lbs / 5.5 kg
Electrical	120V, 50-60 Hz

ORDERING INFORMATION

Catalog No.	Description
C0100-A	Z 100 A Clinical Laboratory Centrifuge with 6 x 15ml rotor
C0200-17A	Combination Adapter Set for 5ml, 7ml, and 10ml tubes, pack of 6



24/7 Customer Service
Phone: (424) 400-2340
Email: info@socalbiomedical.com
www.socalbiomedical.com