



SkySpin Technology

SkySpin CM-6MT - Designed for Versatility

With Multiple Rotors and Adapters

CM-6MT is capable of spinning almost everything in your lab



Compact and Versatile CM-6MT, with its exceptionally slim design, is one of the most compact and versatile swing-out centrifuges on the market.

Most Advanced, Service-Free Drive Powered by the industry's most advanced, service-free drive and used in labs worldwide, the CM-6MT has proven to be an incredibly versatile, long-lasting, and extremely quiet centrifuge that would make an outstanding addition to your laboratory.

Rapid Acceleration & Braking System The CM-6MT has a rotational speed of up to 3500 RPM, and comes with highly responsive rapid acceleration/braking system.

Specifications:	Units	CM6MT
Rotor rotation speed	RPM	100-3500
Maximum force	RCF	2300
Timer	MIN	0.1-99
Working temp range	C	4-40
Relative humidity at 20°C	%	80
Dimensions	mm	430x400x185
Voltage and Frequency	V,Hz	100-240V; 50/60Hz

CM-6MT-8

General centrifuge CM-6MT (rotor not included)

CE

Additional Features:

- Digital Control Panel
- RPM/RCF Indicator
- Sound Notifications
- Digital Rotor Imbalance Sensor
- Cold Room & Refrigerator Safe
- Transperent Lid & Auto Lock



CM-6MT-3 12 Tubes
Maximum volume of applied test tubes: 15 ml.

Maximum sizes of the applied test tubes: (D x L) 16.8 x 140mm.
Top speed: 3500 rpm.



CM-6MT-5 24 Tubes
Maximum volume of applied test tubes: 15 ml.

Maximum sizes of the applied test tubes: (D x L) 16.8 x 115 mm -140mm
Top speed: 3500 rpm.



CM-6MT-4 4 Tubes
Maximum volume of applied test tubes: 50 ml.

Maximum sizes of the applied test tubes: (D x L) 30 x 135 mm.
Top speed: 3500 rpm.



CM-6MT-6 2 Plates
Size of applied plates: (W x L x H) 86 x 128 x 20 mm.

Top speed: 2500 rpm.



CM-6MT-7
5ml Test Tube Adapter.
Compatible with CM-6MT-3 and CM-6MT-5



SkySpin Technology

Clinical SkySpin CM-6MTF - Medical Centrifuge

With Multiple Rotors and Adapters
CM-6MTF is capable of spinning almost everything in your lab



Compact and Versatile CM-6MTF, with its exceptionally slim design, is one of the most compact and versatile FDA registered swing-out centrifuges on the market.

Most Advanced, Service-Free Drive Powered by the industry's most advanced, service-free drive and used in labs worldwide, the CM-6MTF has proven to be an incredibly versatile, long-lasting, and extremely quiet centrifuge that would make an outstanding addition to your laboratory.

Rapid Acceleration & Braking System The CM-6MTF has a rotational speed of up to 3500 RPM, and comes with highly responsive rapid acceleration/braking system.

Specifications:	Units	CM6MTF
Rotor rotation speed	RPM	100-3500
Maximum force	RCF	2300
Timer	MIN	0.1-99
Working temp range	C	4-40
Relative humidity at 20°C	%	80
Dimensions	mm	430x400x185
Voltage and Frequency	V,Hz	100-240V; 50/60Hz

CM-6MT-8 Unit comes with 3 Year Warranty

CE

Additional Features:

- Digital Control Panel
- RPM/RCF Indicator
- Sound Notifications
- Digital Rotor Imbalance Sensor
- Cold Room & Refrigerator Safe
- Transparent Lid & Auto Lock



CM-6MT-3 12 Tubes
Maximum volume of applied test tubes: 15 ml.
Maximum sizes of the applied test tubes: (D x L) 16.8 x 140mm.
Top speed: 3500 rpm.

CM-6MT-5 24 Tubes
Maximum volume of applied test tubes: 15 ml.
Maximum sizes of the applied test tubes: (D x L) 16.8 x 115 mm -140mm

CM-6MT-4 4 Tubes
Maximum volume of applied test tubes: 50 ml.
Maximum sizes of the applied test tubes: (D x L) 30 x 135 mm.
Top speed: 3500 rpm.

CM-6MT-6 2 Plates
Size of applied plates: (W x L x H) 86 x 128 x 20 mm.
Top speed: 2500 rpm.

CM-6MT-7
5ml Test Tube Adapter.
Compatible with CM-6MT-3 and CM-6MT-5