



## 51007 White Polyethylene Core DX Multi Tasking Storage Cart with Clear Door: 23 x 35 x 19" WHD



### Part Number: 51007

Product Dimensions: 23 x 35 x 19 Inches (width x height x depth)

Product Weight: 61 Lbs

Est. Shipping Dimensions: 27 x 39 x 23 Inches

Est. Shipping Weight: 71 lbs

Country of Origin: United States

Suggested Retail Price

**\$579.00**

## Features

- Lab cart with 3 full-extension drawers, 6 side shelves and 1 bulk storage area with clear locking door
- Made of polyethylene, ABS, and clear PETG construction
- Built in push-pull handle in top, 2 keys
- Ideal for organizing, storing and moving any type of supplies for laboratories, hospitals, clinics, industrial labs and others who want an economical, locking cart with clean professional appearance
- 23 x 35 x 19 inches (width, height depth); 61 pounds

## Description

TrippNT's White Polyethylene Core DX Multi Tasking Storage Cart with Clear Door offers an efficient way to sort, store, organize, lock and transport a wide variety of supplies with drawers, shelves, built in glove & tissue dispensers and bulk storage features. Constructed from tough, non-reactive white high density polyethylene, ABS and PETG this lab cart glides smoothly across most surfaces on its 4 heavy duty 3 inch non-marking

rubber casters. The small footprint allows it to be stored securely and unobtrusively out of the way. It features 3 full extension drawers and 6 side compartments, 2 with dispensing slots. A large bulk storage area with a locking sliding, clear PETG door secures supplies. 23 x 35 x 19 inches (width, height depth); 61 pounds