



TrippNT 51142 White Polyethylene Core SP Space Saving Lab Cart with Clear Door: 15 x 35 x 19" WHD



Part Number: 51142

Product Dimensions: 15 x 35 x 19 Inches (width x height x depth)

Product Weight: 52 Lbs

Est. Shipping Dimensions: 19 x 39 x 23 Inches

Est. Shipping Weight: 62 lbs

Country of Origin: United States

Suggested Retail Price

\$559.00

Features

- Compact locking lab cart with 3 full-extension drawers and bulk storage, and clear acrylic door in a small 2 square foot area
- High density polyethylene construction with clear PETG door; 150 pound weight capacity
- 3-Inch non-marking industrial casters; 2 with brakes; built in handle; integrated security lock
- Ideal for mobile storage for any type of supplies work, home, labs, hospitals, clinics or industrial settings
- 15 x 35 x 19 inches (width height depth); 52 pounds

Description

TrippNT's White Polyethylene Core SP Space Saving Lab Cart with Clear Door provides a large amount of extra storage and workspace in a small footprint. This all white cart carries almost three times its weight and the innovative sliding door provides access to your supplies in cramped rooms or tight hallways. The door locks in either the up or down position with an embedded security lock and comes with two keys. Keep supplies

and other equipment in one of the three full-extension drawers or bulk storage area. Built in ergonomic handles make the Core SP easy to maneuver. 15 x 35 x 19 inches (width height depth); 52 pounds