











































## Dispersing Shaft Set Matching &amp; Application Table

Model/ Cat . No.:	Photo Stator	Photo Rotor	Shaft Set Photo	Application Field/Industry	Medium Function	Treatable Volume (H2O)/ml	Tool length (mm)	Min/max emersion depth (mm)	Max viscosity (mPa)	Ultimate suspension fitness ( $\mu$ m)	Ultimate Emulsion fitness ( $\mu$ m)
<b>SET 1</b> <b>H5S2CSR2</b>	 20mm	 15mm		Ceramic Pharmaceutical Petrochemistry Sewage Control	SOLID / LIQUID	10-5000	220	40/175	10,000	10-50	1-10
<b>SET 2</b> <b>H5S2CCR2</b>	 20mm	 15mm		Paper&Tissue Food Tobacco Medicine Sewage Control	STRINGY/ FIBROUS	10-5000	220	40/175	30,000	10-50	1-10
<b>SET 3</b> <b>H5S2CMR2</b>	 20mm	 15mm		Ceramic Paint	SOLID / LIQUID	10-5000	220	40/175	10,000	10-50	1-10
<b>SET 4</b> <b>H5S2FER2</b>	 20mm	 15mm		Paper & Tissue Pharmaceutical Sewage Control Paint	SUSPENSION/ EMULSIONS	10-5000	220	40/175	10,000	10-50	1-10
<b>SET 5</b> <b>H5S2FCR2</b>	 20mm	 15mm		Biotechnology Paper & Tissue Food Tobacco Medicine Sewage Control	STRINGY/ FIBROUS	10-5000	220	40/175	10,000	10-150	1-10

## Dispersing Shaft Set Matching &amp; Application Table

Model/ Cat . No.:	Photo Stator	Photo Rotor	Shaft Set Photo	Application Field/Industry	Medium Function	Treatable Volume (H2O)/ml	Tool length (mm)	Min/max emersion depth (mm)	Max viscosity (mPa)	Ultimate suspension fitness ( $\mu$ m)	Ultimate Emulsion fitness ( $\mu$ m)
<b>SET 6</b> <b>H5S2FMR2</b>	 20mm	 15mm		Ceramic Paper & Tissue Food Cosmetics Petrochemistry Paint	SOLID / LIQUID	10-5000	220	40/175	10,000	10-150	1-10
<b>SET 7</b> <b>H5S3CSR3</b>	 30mm	 23mm		Paper & Tissue Food Pharmaceutical Medicine Sewage Control	SOLID / LIQUID	100-8000	220	40/175	10,000	5-25	1-5
<b>SET 8</b> <b>H5S3CCR3</b>	 30mm	 23mm		Paper & Tissue Food Pharmaceutical Medicine Sewage Control	STRINGY/ FIBROUS	100-8000	220	40/175	30,000	5-25	1-5
<b>SET 9</b> <b>H5S3CMR3</b>	 30mm	 23mm		Ceramic Paint	SOLID / LIQUID	100-8000	220	40/175	10,000	5-25	1-5

## Dispersing Shaft Set Matching &amp; Application Table

Model/ Cat . No.:	Photo Stator	Photo Rotor	Shaft Set Photo	Application Field/Industry	Medium Function	Treatable Volume (H2O)/ml	Tool length (mm)	Min/max emersion depth (mm)	Max viscosity (mPa)	Ultimate suspension fitness ( $\mu$ m)	Ultimate Emulsion fitness ( $\mu$ m)
<b>SET 10</b> <b>H5S3CMR2</b>	 30mm	 15mm		Ceramic Food Sewage Control	MIXING / STIRRING	250-20,000	220	40/175	10,000	HIGH SPEED MIXER	
<b>SET 11</b> <b>H5S3FER3</b>	 30mm	 23mm		Paper & Tissue Pharmaceutical Sewage Control Paint	SUSPENSION/ EMULSIONS	100-8000	220	40/175	10,000	5-25	1-5
<b>SET 12</b> <b>H5S3FSR3</b>	 30mm	 23mm		Paper & Tissue Pharmaceutical Sewage Control Paint	SOLID / LIQUID	100-8000	220	40/175	10,000	5-25	1-5
<b>SET 13</b> <b>H5S3FMR3</b>	 30mm	 23mm		Ceramic Food Cosmetics Pharmaceutical Tobacco	SOLID / LIQUID	100-8000	220	40/175	10,000	5-25	1-5
<b>SET 14</b> <b>H5S4CMR3</b>	 40mm	 23mm		Ceramic Food Sewage Control	MIXING / STIRRING	1000-40,000	220	40/175	10,000	HIGH SPEED MIXER	

### Dispersing Shaft Set Matching & Application Table

Model/ Cat . No.:	Photo Stator	Photo Rotor	Shaft Set Photo	Application Field/Industry	Medium Function	Treatable Volume (H2O)/ml	Tool length (mm)	Min/max emersion depth (mm)	Max viscosity (mPa)	Ultimate suspension fitness ( $\mu$ m)	Ultimate Emulsion fitness ( $\mu$ m)
----------------------	-----------------	----------------	--------------------	-------------------------------	--------------------	---------------------------------	------------------------	--------------------------------------	---------------------------	---	---

**SET 15**

**H5S05CSR04**



Stator : 5mm

Rotor : 4mm

Micro-  
Biotechnology  
Medicine  
Pharmaceutical

SOLID /  
LIQUID

0.2-50

10/60

128

10,000

10-50

1-10