



Cat. No.	Descriptions	HS5-Pro	HS5-S	HS4-Pro	HS4T-Pro	HS3-E	HS2-S
① HS3BCP	Blue carrying plate, use with color quarter pies			■	■	■	■
② HS3BFR	Blue Fixed ring, use with color quarter pies			■	■	■	■
③ HS3RQ	Red quarter pie, 11 holes, 4 mL reaction vessel, Ø15.2mm, 20mm depth	■	■	■	■	■	■
④ HS3PQ	Purple quarter pie, 4 holes, 20 mL reaction vessel, Ø28mm, 24mm depth	■	■	■	■	■	■
⑤ HS3BQ	Blue quarter pie, 4 holes, 30mL reaction vessel, Ø28mm, 30mm depth	■	■	■	■	■	■
⑥ HS2BKQ	Black quarter pie, 4 holes, 40mL reaction vessel, Ø28mm, 43mm depth	■	■	■	■	■	■
⑦ HS2GQ	Green quarter pie, 6 holes, 8mL reaction vessel, Ø17.8mm, 26mm depth	■	■	■	■	■	■
⑧ HS2GDQ	Golden quarter pie, 4 holes, 16mL reaction vessel, Ø21.6mm, 31.7mm depth	■	■	■	■	■	■
⑨ HS2B50	Reaction block for 50mL round bottom flask (one flask capacity)			■	■	■	■
⑩ HS2B100	Reaction block for 100mL round bottom flask (one flask capacity)			■	■	■	■
⑪ HS2B250	Reaction block for 250mL round bottom flask (one flask capacity)			■	■	■	■
⑫ HS2B500	Reaction block for 500mL round bottom flask (one flask capacity)			■	■	■	■
⑬ HS3PC	F101 protective cover, temperature resistance of 135°C			■	■	■	
⑭ HSPT1A	PT1000-A, length of 230mm	■	■	■	■		■
⑮ HSPT1B	PT1000-B, temperature sensor with glass coated, length of 230mm	■	■	■	■		■
⑯ HSPT1SA	Support clamp of PT1000	■	■	■	■		
⑰ HSPT1SB	Support clamp of PT1000						■