



## Tube Mixers

Series	TR6		RT1		RT2	
Style	Roller		Disk		Rotisserie	
Model	TR6-E	TR6-Pro	RT1-E	RT1-Pro	RT2-E	RT2-Pro
Image						
Speed Range	10-70 rpm	10-70 rpm	10-80 rpm	10-70 rpm	10-80 rpm	10-70 rpm
Speed Display	Analog	LCD	Analog	LCD	Analog	LCD
Timer Display	Analog	LCD	Analog	LCD	Analog	LCD
Timer	n/a	1-1199 min	n/a	1-1199 min	n/a	1-1199 min
Operation Mode	Continuous	Continuous / timer	Continuous	Continuous / timer	Continuous	Continuous / timer
Tilt Angle	n/a	n/a	90°	Adjustable, 0-90°	n/a	n/a
Rocking Motion	Rocking and rolling	Rocking and rolling	n/a	n/a	n/a	n/a
Amplitude	24mm	24mm	n/a	n/a	n/a	n/a
Number of Rollers	6	6	n/a	n/a	n/a	n/a
Roller Length	13.5 in	11 in	n/a	n/a	n/a	n/a
Max Tube Diameter	3 in	3 in	n/a	n/a	n/a	n/a
Tube Capacity	n/a	n/a	1.5ml x 48 / 15ml x 16 / 50ml x 12	1.5mL x 60 / 15mL x 16 / 50mL x 8	1.5mL x 12 / 15mL x 12 / 50mL x 12	1.5mL x 48 / 15mL x 24 / 50mL x 24
Included Accessory	n/a	n/a	Disk with 12ea: 1.5/15/50mL clamps	8 x 50mL Disk Accessory	Rotisserie with 12ea: 1.5/15/50mL clamps	16 x 50mL Rotisserie Accessory
Motor Type	DC Brushless	DC Brushless	DC Brushless	DC Brushless	DC Brushless	DC Brushless
Overall Dimensions [L x W x H]	21 x 11 x 5 in	17.5 x 11 x 5 in	9 x 12 x 12 in	10 x 8 x 10 in	21 x 6 x 8 in	21 x 9 x 11 in
Weight (lbs)	12 lbs	9 lbs	6 lbs	10 lbs	6 lbs	10 lbs
Permissible ambient temp and humidity	5~40°C, 80%	5~40°C, 80%	5~40°C, 80%	5~40°C, 80%	5~40°C, 80%	5~40°C, 80%
Protection Class	IP21	IP21	IP21	IP21	IP21	IP21
Voltage	100-240V, 50/60Hz	100-240V, 50/60Hz	100-240V, 50/60Hz	100-240V, 50/60Hz	100-240V, 50/60Hz	100-240V, 50/60Hz
Power	25W	30W	20W	40W	20W	40W
Warranty	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years
Price	\$220	\$506	\$396	\$639	\$396	\$748