


Material : Stainless Steel


Upright Freezer Drawer Racks

For 96 Deep-Well Microtiter Plates



RListed racks are for standard 96 Deep-Well Microtiter Plates and are constructed of corrosion resistant stainless steel, including handles and identification card holders.

Each rack compartment holds one 96 Deep-Well Microtiter Plate. Drawer style rack design allows easy access to fewer plates at a time thus protecting other samples in the rack from warming up. Security lock devices are also available for each rack. Using the internal dimensions obtained for your freezer shelf space (see page 1, fig 1) determine the height (H), length (L) and width (W) of each rack you select for your upright freezer. Racks fit all standard types of upright freezers.



Rack w/DD Lock Device Option

Rack w/96 Deep-Well Plates

Part Number	Description	Dimension L x H x W	Capacity	Lock Device Option*
UFDDP-33	Rack w/ 3 drawers (3 compartments per drawer) for 96 Deep-Well Microtiter Plates	10 ¹⁵ / ₁₆ x 6 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 277 x 154 x 138mm	9	DD-05
UFDDP-34	Rack w/ 4 drawers (3 compartments per drawer) for 96 Deep-Well Microtiter Plates	10 ¹⁵ / ₁₆ x 8 ¹ / ₈ x 5 ⁷ / ₁₆ inch 277 x 205 x 138mm	12	DD-07
UFDDP-35	Rack w/ 5 drawers (3 compartments per drawer) for 96 Deep-Well Microtiter Plates	10 ¹⁵ / ₁₆ x 10 ¹ / ₈ x 5 ⁷ / ₁₆ inch 277 x 256 x 138mm	15	DD-10
UFDDP-43	Rack w/ 3 drawers (4 compartments per drawer) for 96 Deep-Well Microtiter Plates	14 ¹ / ₂ x 6 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 367 x 154 x 138mm	12	DD-05
UFDDP-44	Rack w/ 4 drawers (4 compartments per drawer) for 96 Deep-Well Microtiter Plates	14 ¹ / ₂ x 8 ¹ / ₈ x 5 ⁷ / ₁₆ inch 367 x 205 x 138mm	16	DD-07
UFDDP-45	Rack w/ 5 drawers (4 compartments per drawer) for 96 Deep-Well Microtiter Plates	14 ¹ / ₂ x 10 ¹ / ₈ x 5 ⁷ / ₁₆ inch 367 x 256 x 138mm	20	DD-10
UFDDP-53	Rack w/ 3 drawers (5 compartments per drawer) for 96 Deep-Well Microtiter Plates	18 x 6 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 457 x 154 x 138mm	15	DD-05
UFDDP-54	Rack w/ 4 drawers (5 compartments per drawer) for 96 Deep-Well Microtiter Plates	18 x 8 ¹ / ₈ x 5 ⁷ / ₁₆ inch 457 x 205 x 138mm	20	DD-07
UFDDP-55	Rack w/ 5 drawers (5 compartments per drawer) for 96 Deep-Well Microtiter Plates	18 x 10 ¹ / ₈ x 5 ⁷ / ₁₆ inch 457 x 256 x 138mm	25	DD-10
UFDDP-63	Rack w/ 3 drawers (6 compartments per drawer) for 96 Deep-Well Microtiter Plates	21 ¹ / ₂ x 6 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 547 x 154 x 138mm	18	DD-05
UFDDP-64	Rack w/ 4 drawers (6 compartments per drawer) for 96 Deep-Well Microtiter Plates	21 ¹ / ₂ x 8 ¹ / ₈ x 5 ⁷ / ₁₆ inch 547 x 205 x 138mm	24	DD-07
UFDDP-65	Rack w/ 5 drawers (6 compartments per drawer) for 96 Deep-Well Microtiter Plates	21 ¹ / ₂ x 10 ¹ / ₈ x 5 ⁷ / ₁₆ inch 547 x 256 x 138mm	30	DD-10
UFDDP-73	Rack w/ 3 drawers (7 compartments per drawer) for 96 Deep-Well Microtiter Plates	25 x 6 ¹ / ₁₆ x 5 ⁷ / ₁₆ inch 635 x 154 x 138mm	21	DD-05
UFDDP-74	Rack w/ 4 drawers (7 compartments per drawer) for 96 Deep-Well Microtiter Plates	25 x 8 ¹ / ₈ x 5 ⁷ / ₁₆ inch 635 x 205 x 138mm	28	DD-07
UFDDP-75	Rack w/ 5 drawers (7 compartments per drawer) for 96 Deep-Well Microtiter Plates	25 x 10 ¹ / ₈ x 5 ⁷ / ₁₆ inch 635 x 256 x 138mm	35	DD-10

*Lock device is ordered separately.