



**Material : Stainless Steel**

# Upright Freezer Racks

## For 96 Deep-Well Microtiter Plates







Racks are for 96 deep-well microtiter plates. High quality Upright Freezer Racks, including handles, are constructed of corrosion resistant stainless steel and offer convenient storage of samples in an organized format. Each rack compartment holds a single 96-well deep well microtiter plate. Security lock devices are also available for each rack. Using the internal dimensions obtained for your freezer shelf space (see page 1, fig 1) determine the height (H), length (L) and width (W) of each rack you select for your upright freezer. Racks come with locking rods to secure microtiter plates in place.



Rack w/ Locking Rod



Rack w/DP Lock Device Option

Part Number	Dimension L x H x W	Capacity (Plate)	Lock Device Option*
UFDP-55L	18 <sup>3</sup> / <sub>8</sub> x 9 <sup>3</sup> / <sub>16</sub> x 5 <sup>1</sup> / <sub>2</sub> inch 466 x 233 x 140 mm	5 x 5 =25 	DP-1809
UFDP-56L	18 <sup>3</sup> / <sub>8</sub> x 11 x 5 <sup>1</sup> / <sub>2</sub> inch 466 x 280 x 140 mm	5 x 6 =30 	DP-1911
UFDP-65L	22 x 9 <sup>3</sup> / <sub>16</sub> x 5 <sup>1</sup> / <sub>2</sub> inch 559 x 233 x 140 mm	6 x 5 =30 	DP-2209
UFDP-66L	22 x 11 x 5 <sup>1</sup> / <sub>2</sub> inch 559 x 280 x 140 mm	6 x 6 =36 	DP-2211

\*Lock device is ordered separately.