

Material : Stainless Steel

Upright Freezer Drawer Racks

For Centrifuge Tubes

Racks are specifically designed to provide an ideal alternative for storing large quantities of samples in 15ml and 50ml centrifuge tubes. High quality Upright Freezer Drawer Racks are constructed of corrosion resistant stainless steel, including handles and


identification card holders. Drawer Racks facilitate both access to and storage of samples without having to remove the whole rack from the freezer. Security lock devices are also available for each rack. Using the internal dimensions obtained for your freezer shelf space (see page 1, fig 1) determine the height (H), length (L) and width (W) of each rack you select for your Upright Freezer.





Rack w/50ml Tube



Rack w/DD Lock Device Option



Rack w/15ml Tube

Content to Store	Rack Length	Rack Height	Number of Drawers	Part Number	Description	Dimension L x H x W	Capacity	Lock Device Option*
15ml 	16 1/2 inch 419 mm	4 15/16 inch 126 mm	1	UFD-T15-1	Rack with Drawer for 15 ml Centrifuge Tubes	16 1/2 x 4 15/16 x 5 1/2 inch 419 x 126 x 140 mm	60 Tubes (15 ml)	DD-04
	22 inch 559 mm			UFD-LT15-1		22 x 4 15/16 x 5 1/2 inch 559 x 126 x 140 mm	80 Tubes (15 ml)	
	26 3/4 inch 679 mm			UFD-XLT15-1		26 3/4 x 4 15/16 x 5 1/2 inch 679 x 126 x 140 mm	104 Tubes (15 ml)	
	16 1/2 inch 419 mm	9 15/16 inch 252 mm	2	UFD-T15-2		16 1/2 x 9 15/16 x 5 1/2 inch 419 x 252 x 140 mm	120 Tubes (15 ml)	DD-09
	22 inch 559 mm			UFD-LT15-2		22 x 9 15/16 x 5 1/2 inch 559 x 252 x 140 mm	160 Tubes (15 ml)	
	26 3/4 inch 679 mm			UFD-XLT15-2		26 3/4 x 9 15/16 x 5 1/2 inch 679 x 252 x 140 mm	208 Tubes (15 ml)	
50ml 	16 1/2 inch 419 mm	4 15/16 inch 126 mm	1	UFD-T50-1	Rack with Drawer for 50 ml Centrifuge Tubes	16 1/2 x 4 15/16 x 5 1/2 inch 419 x 126 x 140 mm	30 Tubes (50 ml)	DD-04
	22 inch 559 mm			UFD-LT50-1		22 x 4 15/16 x 5 1/2 inch 559 x 126 x 140 mm	39 Tubes (50 ml)	
	26 3/4 inch 679 mm			UFD-XLT50-1		26 3/4 x 4 15/16 x 5 1/2 inch 679 x 126 x 140 mm	51 Tubes (50 ml)	
	16 1/2 inch 419 mm	9 15/16 inch 252 mm	2	UFD-T50-2		16 1/2 x 9 15/16 x 5 1/2 inch 419 x 252 x 140 mm	60 Tubes (50 ml)	DD-09
	22 inch 559 mm			UFD-LT50-2		22 x 9 15/16 x 5 1/2 inch 559 x 252 x 140 mm	78 Tubes (50 ml)	
	26 3/4 inch 679 mm			UFD-XLT50-2		26 3/4 x 9 15/16 x 5 1/2 inch 679 x 252 x 140 mm	102 Tubes (50 ml)	

*Lock device is ordered separately.