



Material : Stainless Steel

Upright Freezer Racks


For 100-Cell Hinged Top Plastic Storage Boxes




Listed racks are for 100-Cell Hinged Top Plastic Storage Boxes. High quality Upright Freezer Racks, including handles, are constructed of corrosion resistant stainless steel and offer convenient storage of samples in an organized format. Each rack compartment holds one 100-Cell Hinged Top Plastic Storage Boxes. Security lock devices are also available for each rack. Using the internal dimensions obtained for your freezer shelf space (see page 1, fig 1) determine the height (H), length (L) and width (W) of each rack you select for your upright freezer.









Rack Only



Rack w/100-Cell Hinged Plastic Storage Box



Rack w/DR Lock Device Option

Part Number	Description	Dimension L x H x W	Capacity (Box)	Lock Device Option**
UFHT-23	Rack Only	12 1/2 x 6 15/16 x 5 7/8 inch 318 x 176 x 149 mm	2 x 3 =6 	DP-1307
UFHT-23-PX*	Rack w/100-Cell Hinged Plastic Storage Box			
UFHT-24	Rack Only	12 1/2 x 9 1/4 x 5 7/8 inch 318 x 235 x 149 mm	2 x 4 =8 	DP-1309
UFHT-24-PX*	Rack w/100-Cell Hinged Plastic Storage Box			
UFHT-25	Rack Only	12 1/2 x 11 9/16 x 5 7/8 inch 318 x 293 x 149 mm	2 x 5 =10 	DP-1312
UFHT-25-PX*	Rack w/100-Cell Hinged Plastic Storage Box			
UFHT-33	Rack Only	18 3/4 x 6 15/16 x 5 7/8 inch 476 x 176 x 149 mm	3 x 3 =9 	DP-1907
UFHT-33-PX*	Rack w/100-Cell Hinged Plastic Storage Box			
UFHT-34	Rack Only	18 3/4 x 9 1/4 x 5 7/8 inch 476 x 235 x 149 mm	3 x 4 =12 	DP-1909
UFHT-34-PX*	Rack w/100-Cell Hinged Plastic Storage Box			
UFHT-35	Rack Only	18 3/4 x 11 9/16 x 5 7/8 inch 476 x 293 x 149 mm	3 x 5 =15 	DP-1912
UFHT-35-PX*	Rack w/100-Cell Hinged Plastic Storage Box			

* To order replace , X with colour of box : O-orange, P-pink, Y-yellow, B-blue, G-green or N-natural.

**Lock device is ordered separately.