



Material : Stainless Steel

## Upright Freezer Racks

### For 15ml & 50ml Centrifuge Tubes Boxes



Racks are for standard cardboard 15ml or 50ml centrifuge tube boxes. High quality Upright Freezer Racks, including handles, are constructed of corrosion resistant stainless steel and offer convenient storage of samples in centrifuge tubes in an organized format. Security lock devices are also available for each rack. Using the internal dimensions obtained for your freezer shelf space (see page 1, fig 1) determine the height (H), length (L) and width (W) of each rack you select for your upright freezer.



Chest Freezer Rack w/Cardboard Box



w/DR Lock Device Option

Part Number	Description	Dimension L x H x W	Capacity (Box)	Lock Device Option*
UFLB-22	Rack Only	12 <sup>3</sup> / <sub>16</sub> x 10 <sup>3</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>8</sub> inch 310 x 259 x 149 mm	2 x 2 =4	DP-1210
UFLB-22-X*	Rack w/box and divider			
UFLB-32	Rack Only	18 <sup>1</sup> / <sub>4</sub> x 10 <sup>3</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>8</sub> inch 464 x 259 x 149 mm	3 x 2 =6	DP-1810
UFLB-32-X*	Rack w/box and divider			
UFLB-42	Rack Only	24 <sup>3</sup> / <sub>16</sub> x 10 <sup>3</sup> / <sub>16</sub> x 5 <sup>7</sup> / <sub>8</sub> inch 617 x 259 x 149 mm	4 x 2 =8	DP-2410
UFLB-42-X*	Rack w/box and divider			

\* To order , replace X with :16 for 50ml tube box or 36 for 15ml tube box.

\*\* Lock device ordered separately.